

UV PHOTOCHEMICAL OXIDATION AND EXTRACTION OF MARINE DISSOLVED ORGANIC CARBON AT UC IRVINE: STATUS, SURPRISES, AND METHODOLOGICAL RECOMMENDATIONS

Brett D Walker^{1,2*}  • Steven R Beaupré³ • Sheila Griffin² • Ellen R M Druffel²

¹Department of Earth and Environmental Sciences, University of Ottawa, Ottawa, ON K1N 6N5, Canada

²Department of Earth System Science, University of California Irvine, CA 92697-3100, USA

³School of Marine and Atmospheric Sciences, Stony Brook University, Stony Brook, NY 11794-5000, USA

ABSTRACT. The first ultraviolet photochemical oxidation (UVox) extraction method for marine dissolved organic carbon (DOC) as CO₂ gas was established by Armstrong and co-workers in 1966. Subsequent refinement of the UVox technique has co-evolved with the need for high-precision isotopic ($\Delta^{14}\text{C}$, $\delta^{13}\text{C}$) analysis and smaller sample size requirements for accelerator mass spectrometry radiocarbon (AMS ¹⁴C) measurements. The UVox line at UC Irvine was established in 2004 and the system reaction kinetics and efficiency for isolating seawater DOC rigorously tested for quantitative isolation of ~1 mg C for AMS ¹⁴C measurements. Since then, improvements have been made to sampling, storage, and UVox methods to increase overall efficiency. We discuss our progress, and key UVox system parameters for optimizing precision, accuracy, and efficiency, including (1) ocean to reactor: filtration, storage and preparation of DOC samples, (2) cryogenic trap design, efficiency and quantification of CO₂ break through, and (3) use of isotopic standards, blanks and small sample graphitization techniques for the correction of DOC concentrations and Fm values with propagated uncertainties. New DOC UVox systems are in use at many institutions. However, rigorous assessment of quantitative UVox DOC yields and blank contributions, DOC concentrations and carbon isotopic values need to be made. We highlight the need for a community-wide inter-comparison study.

KEYWORDS: dissolved organic carbon, photochemical oxidation, radiocarbon dating, seawater, stable isotope, UV.

INTRODUCTION

Marine dissolved organic carbon (DOC) is the oceans' largest organic carbon reservoir. Ultraviolet photochemical oxidation (UVox) has been used for decades as a method to measure DOC concentrations (noted herein as [DOC]), $\delta^{13}\text{C}$, and $\Delta^{14}\text{C}$ values (Armstrong et al. 1966; Williams 1968; Williams et al. 1969). Renewed interest in DOC has resulted in the development of several custom UVox systems, with commonly sought after methodological improvements that include higher sample throughput (>1× sample per day), smaller C blanks, and smaller sample volumes (Beaupré et al. 2007; Xue et al. 2015). DOC $\Delta^{14}\text{C}$ values measured on systems that reproducibly achieve blank-corrected yields of standards within uncertainty of 100%, $\delta^{13}\text{C}$ and $\Delta^{14}\text{C}$ values within uncertainty of consensus values, and methodological precisions less than natural DOC variability are required to reliably interpret marine DOC biogeochemistry. Therefore, careful, continual examination of all procedures from sample collection to isotopic blank corrections are essential for accurate, reproducible, and intercomparable DOC $\Delta^{14}\text{C}$ measurements performed on each UVox system. In order to minimize uncertainties, parameters including UVox reactor kinetics, trapping efficiencies, manometric precisions, and blank corrections should be characterized routinely (Beaupré et al. 2007). Here, we highlight additional parameters to consider when collecting and preparing samples for UVox. We also share several recent “hurdles” we have encountered as important considerations improving UVox efficiency, quantitative sample C recovery, decreasing carbon blanks and minimizing potential isotopic fractionation within UVox systems.

*Corresponding author. Email: brett.walker@uottawa.ca.

METHODS

UV Photochemical Oxidation at UC Irvine

Marine DOC samples are collected by filtering seawater from Niskin bottles through pre-combusted (540°C, 2 hr) GF/F filters (0.7 µm) held in acid cleaned (10% HCl) 70 mm stainless steel filter holders fitted with acid cleaned platinum cured silicone tubing. All GF/F filters are pre-combusted in individual Al foil “envelopes” to minimize C contamination during sample collection. The filtered seawater is collected in pre-combusted (540°C, 2 hr) 1-L amber Boston round bottles. Each bottle is rinsed 3 times with 100 mL of sample prior to filling with ~800 mL of seawater, and then sealed with a Teflon lined cap. An additional layer of chromic acid cleaned PTFE tape is placed between the cap and bottle. Samples are immediately frozen at sea in a 20°C chest freezer, at a 15–20° angle to prevent bottle breakage during freezing.

For detailed discussion of most DOC line run parameters, including a diagram of our vacuum line and reactor systems, we refer the reader to Beaupré et al. (2007). We routinely use the dilution method for UVox isolation of DOC from seawater (Beaupré et al. 2007; Griffin et al. 2010; Druffel et al. 2013; Walker et al. 2016). Here, seawater DOC samples (~800 mL) are diluted to ~1000 mL with 18.2 MΩ cm Milli-Q water (MQ; 0.4–1.0 µM C) prior to UVox. There are several advantages in using the dilution method. First, it effectively reduces the detection limit of a DOC line instrument to very small volumes of seawater (30 mL) (Beaupré et al. 2007). Second, this method saves time by eliminating the need to pre-irradiate MQ water to be used as diluent for samples and standards. It also reduces the total number of standards that are needed for AMS corrections, and thus preparation time from two days to one day (Griffin et al. 2010). Third, dilution allows for routine monitoring of MQ blank C and DOC system line blanks that can be directly used to correct sample concentrations and isotopic values via mass balance. Finally, it prevents evaporative loss of seawater and a buildup of salts on the reactor wall closest to the lamp—preventing a decrease in UVox efficiency.

During sample loading, diluent and sample water volumes are determined non-intrusively by measuring the meniscus height from outside the quartz reactor using a homemade, calibrated cathetometer. The cathetometer is comprised of a 30 cm ruler fixed vertically, and level in both x and y, to a laboratory ring stand. Reactor meniscus heights are read using a 3× power hand lens to an accuracy of 0.5 mm. Reactor-specific volume (mL) and meniscus height (cm) relationships are experimentally determined by measuring known aliquots of MQ water by mass and applying a least squares linear regression analysis. MQ water density is also inferred using a linear density approximation (0.9982 g/mL at 20°C, 0.9957 g/mL at 30°C) after measuring its temperature. The resultant linear regression is then used to determine daily sample and diluent volumes based on their meniscus heights with propagated uncertainties of ±3–5 mL.

After loading, the sample is then acidified with 1 mL 85% phosphoric acid (Fisher HPLC grade). Dissolved inorganic carbon (DIC) is sparged from the sample using ultra high purity (UHP; 5.0) He gas. The mixing ratio of CO₂ venting from the sample is monitored with an infrared gas analyzer (LiCOR, LI-6252) and logged to a lab computer. The LiCOR is calibrated bi-monthly with standard reference gases (0–500 ppm CO₂ in He). Following complete DIC removal, DOC samples are UV oxidized to CO₂ using a 1200W medium-pressure mercury arc lamp for 4 hr. The CO₂ produced photochemically from DOC is then sparged with UHP He, purified cryogenically, and manometrically quantified.

At all stages, UHP He flow rates are quantified with a 50 mL soap-film bubble flow meter (Restek, #RES-20136) at the LiCOR's vent.

Following manometry, equilibrated gas splits of CO₂ are taken for Fm and δ¹³C measurements at UCI's Keck Carbon Cycle AMS (KCCAMS) Laboratory. One split (2–800 μgC) is graphitized via sealed-tube Zn reduction (Xu et al. 2007; Walker and Xu 2019) for accelerator mass spectrometry radiocarbon (AMS ¹⁴C) measurements. The other split is cryogenically transferred to 3 mm diameter, 5.5 cm length, Pyrex tubes for δ¹³C analysis with a Gas Bench II and Finnigan Delta Plus IRMS (Walker et al. 2017). Consensus values for isotopic standards are either known (e.g. from IAEA) or experimentally determined via elemental analysis and IRMS (EA-IRMS; ±0.1‰).

The total reproducibility of these methods for [DOC], Fm, and δ¹³C values are [DOC] ±1.3 μM, Fm ±0.004, and δ¹³C ±0.2‰, respectively (Walker et al. 2016). However, in some cases we have observed Fm of ±0.010 when seawater samples were split in half and analyzed after long storage periods (Druffel et al. 2013). The system “line” blank (i.e. re-irradiated MQ water) ranges from 3 to 12 μgC, but is typically ~4 μgC.

RESULTS AND DISCUSSION

Collection, Storage, and Preparation of Seawater DOC Samples

Extreme care should be taken during DOC sample collection, storage, and laboratory preparation to minimize any possible extraneous C addition. Previous work has shown that acidification and room temperature storage of filtered seawater results in loss of semi-labile (¹⁴C modern) DOC (Walker et al. 2016). While frozen storage does not appear to be affected, it is prone to the unpredictable precipitation of salt crystals that can introduce [DOC] and isotopic measurement errors (Beaupré and Druffel 2009). For example, Druffel et al. (2013) noted changes in DOC recovery (by 1–5 μM), Δ¹⁴C offsets (±10‰, or Fm ±0.0100), and δ¹³C bias (–0.3‰) for duplicate aliquots of seawater from the same sample bottles that were measured on different days. Thus, phase transitions and inconsistent transfer of sample carbon from storage bottles to the UVox reactor may lead to measurement artifacts (Beaupré and Druffel 2009).

In order to minimize errors associated with frozen storage, we originally restricted our measurements to samples that did not develop visible crystals (Beaupré and Druffel 2009). To increase the proportion of useable samples, Walker et al. (2016) quantitatively transferred the contents of each sample bottle into the reactor via rinsing to achieve reproducibility in [DOC], Fm, and δ¹³C of ±1.3 μM, ±0.0094, and ±0.1‰, respectively. Since these crystals are likely precipitated carbonates (Beaupré et al. 2009), we now acidify and shake thawing samples in their storage bottles to re-dissolve crystals before transferring seawater to the reactor (Figure 1). The acidified seawater and any remaining crystals are then transferred quantitatively to the reactor with three or four ~50 mL MQ rinses. These simple procedures have provided reliable storage, compositional fidelity, and decreased variability in [DOC], Δ¹⁴C (or Fm) values (Griffin et al. unpublished data).

Liquid Nitrogen Trap Efficiency and Minimizing Break Through (BT)

The DOC line at UC Irvine utilizes two highly efficient compact Horibe traps to cryogenically trap water (–78°C, dry ice and isopropanol slush) and sample CO₂ (–196°C, liquid nitrogen);




Figure 1 Dissolution of crystals from thawed seawater samples. Pictured is a time series of a thawed, acidified sample with abundant crystals. The sample was vigorously shaken every 1 min to achieve crystal re-dissolution. (A) 1.0 mL of phosphoric acid is added to the 1-L amber Boston round sample bottle, a flashlight behind the bottle helps illustrate the mound of crystals. (B) Half of the original crystals remain after 6 min of repeated shaking. (C) Almost 100% of crystals have re-dissolved by 20 min of shaking.

LN2) (Beaupré et al. 2007). In early tests of the DOC line, He streams with high CO_2 mixing ratios and flow rates were found to allow some sample CO_2 to “break through” (BT) the LN2 Horibe trap and leave the system without being quantified. CO_2 BT was experimentally determined to be minimized through an initial reduction in He flow rate (to below 50 mL min^{-1}) during the sparging of sample DOC. However, data from these experiments were never published and unfortunately, the procedure of reducing the DOC initial sparge rate, was never adopted as standard lab practice. Thus, throughout much of the DOC line’s history (2008–2017), BT was present and some sample C was lost from our DOC samples. An example of a LiCOR trace during a full day experiment is shown in Figure 2 and illustrates this C BT during the DOC sparge step. We have conducted a series of experiments, seeking to quantify and minimize C BT.

Quantifying BT as a Function of [DOC] and Flow Rate

We analyzed 8 OX-I standards to determine the amount of BT as a function of [DOC] in the reactor, and the flow rate of He gas during the DOC sparge. First, BT was measured at two flow rates of He (200 mL min^{-1} and 250 mL min^{-1}) during the DOC sparge, at [DOC] values ranging from 20 to $75 \mu\text{M}$. Results show that BT at both flow rates are equal ($22 \pm 2 \mu\text{gC}$) for samples with [DOC] values of 12 and $15 \mu\text{M}$ (Figure 3A). However, for samples with higher [DOC] values ($>30 \mu\text{M}$), BT is consistently higher for samples that were sparged at 250 mL min^{-1} than those sparged at 200 mL min^{-1} . The average of three BT values for tests run at 200 mL min^{-1} is $27.3 \pm 2.9 \mu\text{gC}$, and that for tests run at 250 mL min^{-1} is $34.9 \pm 0.4 \mu\text{gC}$. We observe no significant differences between Fraction Modern (Fm) values of samples run at the different He flow rates or [DOC] values (Figure 3B). Stable C isotopic ($\delta^{13}\text{C}$) results also show no significant differences in the samples (Figure 3C). In this limited dataset, we do not observe fractionation of carbon isotopes with changing [DOC] or He flow rate.

Following these findings, we have made corrections to our [DOC] values for samples to account for the lost CO_2 (BT) run at UC Irvine prior to these tests; this has resulted in a 3–4 μM increase in our [DOC] values. In addition, we have adopted the practice of DOC sparging at a He flow rate of 45–54 mL min^{-1} . This practice has minimized our DOC sample C BT to $<3 \mu\text{gC}$ for samples and isotopic standards (200–700 μgC), and 0.5 to 1.0 μgC for MQ blanks (10–15 μgC) and line blanks (typically $\sim 4 \mu\text{gC}$). By adopting these




Figure 2 Example of monitoring CO_2 gas leaving the DOC UVox line over an experiment. At the beginning of the day, DIC and residual CO_2 in the DOC line is removed via acidification of the sample and sparging with UHP He (DIC Sparge). When all residual CO_2 is removed, the He flow is stopped, the line isolated from atmosphere and UVox started. Throughout UVox, the system vent tubing slowly “bleeds up” to lab air. This does not affect the DOC sample, because it is isolated. After 4 hr, the UVox is stopped and He gas is used to sparge all resultant CO_2 for isolation in the Horibe trap (DOC Sparge). In this example, a flow rate of 200 mL min^{-1} was used. Several minutes later, we see sample C breaking through the Horibe trap and leaving the system (0–10 ppm). BT lasts between 45 and 70 min, depending on [DOC] and flow rate.

changes (i.e. adding C lost via BT and subtracting the line blank), we routinely achieve DOC standard yields of $100 \pm 2\%$.

Testing of New Stainless Steel Cryogenic Trap Designs

Following the discovery of Horibe trap C BT, we developed and tested several new cryogenic trap designs. One hypothesis for the observed C BT is that as sample CO_2 builds up at the LN2 interface of the glass Horibe trap wall, it may reach a critical mass or geometric shape that permits CO_2 “snow” to blow away, sublime and leave the system with the flowing He gas. One economical solution to this problem was to build a stainless steel trap containing many loops, forcing sample CO_2 to pass through multiple LN2 and room temperature interfaces. Stainless steel has the added benefits of being easy to work with and a much higher thermal conductivity ($16 \text{ W m}^{-1} \text{ K}^{-1}$) than glass ($1 \text{ W m}^{-1} \text{ K}^{-1}$) (Lide 2005), and thus should be more efficient at trapping CO_2 . We have tested three types of stainless steel cryogenic traps for sample C BT.

The first trap was made from 0.25 inch OD 304 stainless steel tubing (3.5 m length, ID=5.3 mm) wrapped to form 7 loops ~6 inches in height. The trap was adapted to the DOC line with Swagelok Ultra-Torr fittings and vacuum tested. Only the bottom third of the trap was placed in LN2, thus forcing sample CO_2 to cross a total of 14 thermal interfaces between room temperature (20°C) and LN2 (-196°C). This trap was tested via acidification of bicarbonate (DOC = $800 \mu\text{M}$) in MQ and He stripping at 200 mL min^{-1} to simulate a DOC sample (as above in the section *Quantifying Break Through (BT) as a Function of [DOC] and Flow Rate*). In this test, we quantified $70 \pm 7 \mu\text{gC}$ BT in the first hour of sparging. A




Figure 3 Mass and C isotopic effects of sample breakthrough (BT) at differing OX-I [DOC] and flow rates. (a) relationship between mass of BT as CO_2 (μgC) at two He flow rates and three [DOC] values. (b) DOC OX-I Fm values, both uncorrected and corrected, showing no effect of sample BT C on corrected ^{14}C measurements. (c) uncorrected DOC OX-I $\delta^{13}\text{C}$ values show no isotopic fractionation with sample C BT. In plots B and C, the green line represents the consensus Fm and $\delta^{13}\text{C}$ value of OX-I. (Please see electronic version for color figures.)

second test (DOC = 800 μM) of this trap with 6-mm Pyrex rods installed, to force sample gas onto the side walls of the trap, yielded a nearly identical amount of C BT (78 ± 7 μgC).

We tested a second trap design made from 0.125 inch OD 304 stainless tubing (3.6 m in length, ID=1.4 mm) with a total of 8 loops and 16 thermal interfaces. In our first bicarbonate




Figure 4 Hybrid stainless steel cryogenic trap design and C BT results. (A) Side-view of the hybrid 0.25 in and 0.125 in OD stainless cryotrap next to a 15 cm ruler. (B) View of hybrid stainless trap in use. The bottom third of the cryotrap is placed in liquid nitrogen (LN2) allowing CO₂ gas to heat up and re-freeze on each pass around the trap. (C) LiCOR trace of CO₂ (ppm; black line) and rate of CO₂ loss (µgC min⁻¹; green line) exiting the system as CO₂ (400 ppm) in He reference gas is flowed at 150 mL min⁻¹. (D) LiCOR trace of CO₂ (ppm; black line) and rate of CO₂ loss (µgC min⁻¹; green line) exiting the system as CO₂ (400 ppm) in He reference gas is flowed at 50 mL min⁻¹. Plots C and D, the blue line represents the baseline CO₂ ppm value achieved with pure UHP He (no CO₂) prior to the beginning of the experiment. Black arrows indicate the moment liquid nitrogen (LN2) was added to the trap, at which point the resultant LiCOR CO₂ ppm values decrease. Grey shaded regions show the 10 min integration periods, from which equilibrated trap BT CO₂ ppm values were determined.

(DOC = 800 µM acidification/sparge test, the trap clogged within the first 5 min of the experiment, likely due to the small tubing diameter, high CO₂ trapping efficiency and the large mass of C being sparged. We observed zero C BT during this experiment.

A third, hybrid trap design, was tested where sparged CO₂ was first passed through the 0.25 in trap, and next via the 0.125 in trap coupled with a Swagelok Ultra-Torr fitting. This trap had a path length of 7.1 m, a volume of 84 mL and a combined total of 15 loops and 30 thermal interfaces (Figure 4A,B). Instead of using bicarbonate acidification, we performed two tests with a standard reference gas (400 ppm CO₂ in UHP He) at 150 mL min⁻¹ and 50 mL min⁻¹. These tests resulted in similar rates of C BT (0.020 µgC min⁻¹) at both flow rates (Figure 4C,D). Once new baseline values were reached for the reference gas flowing through the LN2 trap, we determined a total C BT mass of 2.08±0.21 and 2.30±0.23 µgC at the low and high flow rate condition through integrating 10 minutes of baseline, respectively. Rates of C BT were slightly higher than the glass Horibe trap at 50 mL min⁻¹




Figure 5 Conceptual diagram of carbon mass contributions and aqueous volumes for DOC samples, standards, diluent and line blanks. Various amounts of extraneous C are added for each sample type and diluent volume used. A ubiquitous, volume-less extraneous carbon “line blank” is inherent to all DOC UVox measurements. A seawater sample prepared using the dilution method also contains 100–200 mL of low carbon Milli-Q water as a diluent. A DOC isotopic standard (0.2–0.8 mgC) is accurately weighed and quantitatively transferred to the reactor with diluent. Standards contain far more MQ diluent by both mass (10–15 µgC) and aqueous volume (1000 mL) than samples. A diluent blank (MQ) is comprised only of MQ carbon and the line blank.

(0.001 µgC min⁻¹). Since our modified trap designs did not appear to yield better results, we currently recommend using a glass Horibe trap (design specified in Beaupré et al. 2007), at a DOC sparge flow rate of 50 mL min⁻¹ to minimize C BT. In the future, we intend to test copper tubing traps, which are more affordable, and have a much higher thermal conductance (~400 W m⁻¹ K⁻¹) (Lide 2005).

DOC Line Extraneous C Sources and Use of Standards for Blank Corrections

There are two methods for determination of extraneous C blanks and correction of ¹⁴C data in a DOC UVox system (Figure 5). The first is the *direct method*, where extraneous C blank masses and Fm values are experimentally measured and then these C contributions subtracted from sample measurements via isotopic mass balance. With the direct method, it is imperative that this isotopic mass balance include propagation of known measurement uncertainties for both samples and extraneous C blanks. The second method is the *indirect method* where extraneous C blanks are quantified indirectly by measuring small to ultra-small sizes of both ¹⁴C-free and modern standard reference materials. The ¹⁴C-free standards are used to determine masses of extraneous modern carbon (MC; Fm=1) contributions, and modern standards (OX-1) for extraneous dead carbon (DC; Fm=0) contributions. These DC contributions are also commonly referred to as C backgrounds. The indirect method and mass balance equations are clearly presented in Santos et al. (2007). With the indirect method, ±50% uncertainty for both MC and DC mass contributions are prescribed as a precaution.

Both the *direct* and *indirect* methods have their own advantages and disadvantages for DOC data. For example, the indirect method does not inherently require routine monitoring of MQ diluent or line blank C masses or Fm values (albeit we strongly recommend users do this). The indirect method also does not require ultra-small mass graphitization, which can be difficult for

many labs to achieve and measure on an AMS. However, there are a few disadvantages to the indirect method for DOC line data: (1) it requires analysis of many analytical “process” isotopic standards (ideally $3\times$ for MC and $3\times$ for DC varying in size between 50 and 500 μgC) in order to accurately determine Fm blank contributions, (2) standards are not fully representative of sample matrix (for seawater), and also overestimate MQ C (10–15 μgC) contributions to samples containing <3 μgC (MQ), and (3) this method does not adequately control for the presence of a ubiquitous DOC system line blank.

There are many reasons users should move towards the *direct method* of blank correction, because it: (1) can be more easily implemented given recent advances in ultra-small mass (2–25 μgC), sealed tube Zn graphitization (Walker and Xu 2019), (2) does not overestimate C contributions from MQ diluent, (3) can effectively correct for ubiquitous system line blanks underlying all UVox measurements, (4) yields more representative [DOC] and isotopic values with propagated uncertainties, and (5) can be applied with sample matrix matching (i.e. the addition of pre-combusted salt to seawater standards and blanks), which should control for possible UVox fractionation due to halides in solution. We note that, in the case of UVox systems, the *direct method* requires frequent analysis of both MQ diluent and line blanks to constrain blank Fm variability. This variability can be a function of changes in user operation, system parameters, MQ [DOC], and extraneous organic C sources, such as upstream CO_2 sorbents (e.g. Ascarite II), and the KI trap solution.

Direct measurement of ultra-small DOC line blanks and diluents has historically precluded the implementation of the direct method for UVox systems. At UC Irvine, the vast majority of our samples have been corrected using the *indirect method*. Recently, our group has started implementing *direct method* corrections through quantification of MQ diluent and DOC line blank Fm values using the ultra-small mass sealed tube Zn method (Walker and Xu 2019). In the case of DOC line blanks, the validity of a 0 μgC contribution from re-irradiated MQ (RIMQ) ultimately depends on the handling/storage of this re-irradiated MQ to prevent addition of extraneous C. We measure MQ and RIMQ on sequential UVox days, keeping the RIMQ in the quartz reactor overnight, in order to minimize contamination of the line blank. We routinely measure both modern and ^{14}C -free isotopic standards (OX-I and Glycine, respectively), with a MQ, line blank, or isotopic standard run approximately every 6–8 samples (bi-weekly). We note that salt is currently not added to the diluent used for our isotopic standards. This may affect *indirect-method* blank corrections as standards in a saline matrix may behave differently than in fresh water during UVox would likely have a different extraneous C blank.

Historical Record of DOC UVox Isotopic Standard Fm Values at UC Irvine

Corrected Fm and $\Delta^{14}\text{C}$ values are typically reported for DOC samples. However, to date only a handful of studies have reported percent yields and ^{14}C values of standard reference materials run on DOC UVox systems (Beaupré et al. 2007; Griffin et al. 2010; Beaupré and Druffel 2012; Druffel et al. 2013). Since the establishment of the DOC UVox method at UC Irvine, we have measured over 170 isotopic standards for the correction of [DOC] and Fm values. A summary of corrected and uncorrected Fm values for these isotopic standards is shown in Figure 6. Note that DOC standards run in 2017–2018 were blank corrected using the *direct method*.

Over the past decade, the majority of our measured DOC OX-I standard Fm values are correctable within 2-sigma of the consensus OX-I value (Fm=1.0398). DOC OX-I corrected with the *direct method* (2017–2018) are less variable and more accurate (closer to




Figure 6 Fraction modern (F_m) values of DOC isotopic standards at UC Irvine. (A) DOC OX-I F_m values ($n=109$) vs. date of UVox, (B) DOC OX-I F_m values vs. size, (C) DOC Glycine F_m values vs. date of UVox, (D) DOC Glycine F_m values ($n=67$) vs. size. In all plots, open and closed symbols represent uncorrected and blank corrected F_m data, respectively. Green lines indicate OX-I and Glycine consensus F_m values. Error bars represent 1-sigma propagated errors of blank corrected F_m values.

the consensus OX-I F_m value). This suggests an improved dead C blank correction over the *indirect method* (i.e. less extraneous MQ C is subtracted from samples). While propagated uncertainties increase with decreasing DOC OX-I sample size, F_m values remain correctable within 2-sigma propagated error of the consensus OX-I F_m value (Figure 6B). Overall, our ^{14}C -free DOC Glycine F_m values (2012–present) are more variable (Figure 6C,D), suggesting DOC samples are more susceptible to modern extraneous C contributions. Of the few DOC Glycine samples we have measured in 2017–2018, the *indirect method* also slightly under-corrects for modern C contamination. We typically achieve a combusted Glycine F_m value <0.0020 at UC Irvine, however many corrected DOC Glycine's have F_m values >0.0050 . Aside from larger propagated uncertainties with smaller standard sizes, there does not appear to be a trend between corrected DOC Glycine F_m values and sample size (Figure 6D).

In order to evaluate the difference between *direct* vs. *indirect* blank corrections on DOC F_m values, we use a Bland-Altman (difference) plot on a dataset of $n=34$ DOC samples in which both corrections were applied (Figure 7). The coefficient of determination of least squares regression is weak ($R^2=0.26$), however it is statistically significant, with slope errors excluding zero ($m=0.0012\pm 0.0004$). These results show that modern samples ($F_m=1.0$) are less affected




Figure 7 Bland-Altman difference plot illustrating direct and indirect blank corrected DOC sample Fm data. Green line indicates no difference between the two correction methods. The dashed blue line represents the least squares linear regression of the method difference vs. average Fm value of the two methods. Data are from $n=34$ DOC sample “unknowns” measured at UC Irvine in March 2018.

by the blank correction methods than half-modern ($F_m=0.5$) or dead samples ($F_m=0.0$). This indicates that the *indirect method* is over correcting for MQ C contributions ($F_m=0.4-0.7$). Direct blank corrected Fm values are lower than indirect corrected value by 0.0021 ± 0.0006 ($\Delta^{14} = 2.1 \pm 0.6\text{‰}$). In addition, propagated Fm individual DOC measurement uncertainties with the *direct method* ($F_m \pm 0.0052\text{‰}$) are more representative of our true analytical uncertainty ($F_m \pm 0.0040-0.0100$, or $\Delta^{14} \pm 4-10\text{‰}$), than individual measurement errors determined via the *indirect method* ($F_m \pm 0.0020$).

Historical Record of DOC UVox Isotopic Standard $\delta^{13}C$ Values at UC Irvine

Most UVox studies have focused primarily on [DOC] and ^{14}C measurements. Wide ranges in seawater DOC $\delta^{13}C$ values ($\sim 3\text{‰}$) have been reported from similar depths and locations (e.g. Broek et al. 2017; Zigah et al. 2017). Studies seeking to constrain UVox DOC $\delta^{13}C$ values have received less attention. Here we discuss the historical record of DOC isotopic standard $\delta^{13}C$ values measured at UC Irvine (Figure 8).

Since 2008, our DOC OX-I $\delta^{13}C$ values have ranged from -18.8‰ to -22.0‰ with an average value of -19.5‰ ($\pm 0.5\text{‰}$ standard deviation, $n=95$; Figure 8A). Since 2012, DOC Glycine $\delta^{13}C$ values range from -38.1‰ to -40.8‰ with an average value of -40.1‰ ($\pm 0.5\text{‰}$ standard deviation, $n=46$; Figure 8C). Prior to October 2013, we did not fully equilibrate DOC CO_2 gas splits for $\delta^{13}C$. Instead, sample gas was “bled” off the main sample through a stopcock. We hypothesize that much of the $\delta^{13}C$ variability observed for both OX-I and Glycine (average $\delta^{13}C = -19.6 \pm 0.6\text{‰}$, $n=55$ and $-40.0 \pm 0.9\text{‰}$, $n=9$, before October 2013, respectively) can be attributed to mass-dependent isotopic fractionation induced by this process. Small samples would display a higher degree of mass-dependent fractionation (Figure 8B,D).




Figure 8 Stable isotopic ($\delta^{13}\text{C}$) permil values of DOC isotopic standards at UC Irvine. (A) DOC OX-I $\delta^{13}\text{C}$ values vs. date of UVox, (B) DOC OX-I $\delta^{13}\text{C}$ values vs. samples size, (C) DOC Glycine $\delta^{13}\text{C}$ value vs. date of UVox, (D) DOC Glycine $\delta^{13}\text{C}$ value vs. sample size. In all plots, open symbols represent uncorrected $\delta^{13}\text{C}$ values. Green lines indicate consensus $\delta^{13}\text{C}$ values for OX-I (IAEA reported value) and Glycine (determined at UCI by routine EA-IRMS and closed tube combustion measurements). Error bars represent assigned $\pm 0.2\text{‰}$ errors of UVox $\delta^{13}\text{C}$ values.

After October 2013, we began fully equilibrating CO_2 gas splits by expanding into a larger volume with a 2 min equilibration time prior to “splitting” for ^{13}C ($\sim 12\text{--}35\ \mu\text{gC}$), ($\sim 300\text{--}500\ \mu\text{gC}$) and an archive ($100\ \mu\text{gC}$). With gas equilibration, DOC OX-I and Glycine values are less variable and more accurate, with $\delta^{13}\text{C}$ values 0.1‰ closer to consensus values ($-19.5 \pm 0.2\text{‰}$, $n=40$ and -40.1 ± 0.4 , $n=37$, respectively). While our recent BT tests did not reveal apparent $\delta^{13}\text{C}$ fractionation (see above in section *Quantifying BT as a Function of [DOC] and Flow Rate*), excluding $n=1$ OX-I and Glycine outliers, over the past year our $\delta^{13}\text{C}$ accuracy has overall improved and our variability decreased after controlling for C BT ($-19.3 \pm 0.1\text{‰}$, $n=13$ and $-40.2 \pm 0.1\text{‰}$, $n=7$). Currently, our isotopic standards are within 0.2‰ and 0.4‰ of consensus $\delta^{13}\text{C}$ values for OX-I and Glycine, respectively. October 2013, we have prescribed 1-sigma $\pm 0.2\text{‰}$ DOC $\delta^{13}\text{C}$ errors to help bracket these offsets.

DOC OX-I $\delta^{13}\text{C}$ values are almost always equal to or lower than consensus EA-IRMS values ($-19.1 \pm 0.1\text{‰}$), whereas DOC Glycine's are almost always higher than consensus EA-IRMS values ($-40.6 \pm 0.1\text{‰}$). Possible explanations for these offsets include (1) varying contributions of extraneous C blanks from addition of MQ water as a diluent, (2) intramolecular isotopic heterogeneity of $^{13}\text{C}/^{12}\text{C}$ ratios (e.g. the amino-C vs. carboxyl-C groups

in Glycine have different $\delta^{13}\text{C}$ values), and (3) the UVox process induces C-isotope fractionation differentially for these two compounds. While we have not tested the latter two processes, we believe MQ and extraneous C blank contributions could largely explain our observed offsets (i.e. in both Fm and $\delta^{13}\text{C}$). We are currently working towards examining these effects more closely. For example, we have recently measured $\delta^{13}\text{C}$ of our MQ diluent to be -26‰ (Walker et al. unpublished data). Thus, MQ C contributions to DOC standards would make DOC OX-I $\delta^{13}\text{C}$ more negative and DOC Glycine more positive than their canonical values. The $\delta^{13}\text{C}$ value of OX-I is $\sim 7\text{‰}$ different from that of MQ and Glycine is $\sim 14.6\text{‰}$ different from that of MQ. Given the mass contribution of MQ C to DOC standards and these isotopic endmember offsets, DOC Glycine values should therefore be more affected by blank C contributions than OX-I. This is consistent with our observed offsets of DOC standards from consensus values ($\Delta\delta^{13}\text{C} = -0.2\text{‰}$ and $+0.4\text{‰}$ for OX-I and Glycine, respectively). Small DOC standards should be more affected by MQ C contributions, because they comprise a larger relative percent mass of the DOC standard. However, we do not currently have enough data to make this comparison. Finally, we still working towards determining the best method for DOC $\delta^{13}\text{C}$ blank correction. In the interim, our lab reports “uncorrected” DOC $\delta^{13}\text{C}$ values with prescribed 1-sigma errors of $\pm 0.2\text{‰}$.

UVox Community and Methodological Recommendations

Based on our experience, we have several methodological recommendations for consideration by the broader UVox community. These are especially important for users developing their own homebuilt systems, or using adapted commercial UVox systems (e.g. Ace Glass #7900). First, effectively parameterize the UVox efficiency and kinetics of the system. Particular attention should be paid to the system lamp geometry, power, and photon flux reaching the reactor. In addition, reactor design, sample homogenization (via stirring), cooling and vacuum line design are important considerations. For detailed discussion of these parameters, we refer the reader to Armstrong et al. (1966), Beaupré et al. (2007), and Oppenländer (2007).

For sample preparation and UVox analysis, we recommend care to minimize contamination and maximize procedural efficiency and consistency. Samples should be filtered using pre-combusted GF/F filters ($0.7\ \mu\text{m}$) and stored frozen, as opposed to acidified at room temperature (Walker et al. 2016). When loading samples into the reactor, it is imperative that all crystals (if present) be loaded into the reactor by rinsing with MQ diluent prior to UVox. The presence of crystals precludes the loading of partial samples from a single bottle into the reactor. A low [DOC] diluent should be used (e.g. Milli-Q) and its C content measured regularly for proper correction of [DOC] and isotopic values. Any CO_2 leaving the UVox system, including that which breaks through the LN2 trap during DOC sparge, should be quantified and used to correct reported [DOC] data. This step requires the use of a sensitive CO_2 analyzer (0–400 ppm range).

We recommend routine analysis of diluent blanks, line blanks and isotopic standards, such that quantitative (100%) yields of standards and samples can be ensured and appropriate direct or indirect corrections can be made to DOC $\Delta^{14}\text{C}$ values. At UC Irvine, we typically measure the following standards and blanks interspersed between sets of six DOC samples; a MQ diluent and subsequent DOC line blank (RIMQ) the following day, a large DOC OX-1 (500 μgC), a small DOC Glycine (200 μgC), another set of MQ diluent/line blanks, a small DOC OX-I (200 μgC) and a large DOC Glycine (500 μgC).

For the radiocarbon community, we recommend that UVox users report at a minimum: AMS lab numbers, [DOC], Fm or $\Delta^{14}\text{C}$ of samples, standards (with % yields) and blanks, the year of collection, year of measurement, and IRMS determined $\delta^{13}\text{C}$ values (if AMS ^{13}C values are measured, state this but do not report values). In addition, propagated [DOC] and isotopic uncertainties must be determined and reported. With many new UVox systems coming online, the need for a community-wide inter-comparison study is apparent. UVox measurements should result in *quantitative yields* (i.e. 100%), after blank and BT corrections are applied and errors propagated. A UVox system producing low yields will fractionate DOC based on its inherent photochemical lability. Here, a UVox measurement will not be representative of total DOC, nor directly comparable with measurements from other UVox systems. Publishing standard yields and isotopic data will help users assess the inter-comparability of results—thus helping to shape future interpretations of DOC cycling in the environment.

ACKNOWLEDGMENTS

We acknowledge Noreen Garcia for help with graphitization, and general laboratory assistance. Jennifer Walker, Dachun Zhang and Xiaomei Xu aided with $\delta^{13}\text{C}$ analysis at UC Irvine. Sample Fm values were determined by us at the UC Irvine W. M. Keck Carbon Cycle Accelerator Mass Spectrometry Laboratory. We thank John Southon for his advice and help with AMS analysis. This work was funded by NSF Chemical Oceanography program (OCE-141458941 to E.R.M.D., and OCE-1536597 to S.R.B.) and an American Chemical Society Petroleum Research Fund New Directions (PRF-55430-ND2) Grant (to B.D.W. and E.R.M.D.).

REFERENCES

- Armstrong FA, Williams PM, Strickland JD. 1966. Photo-oxidation of organic matter in sea water by ultra-violet radiation analytical and other applications. *Nature* 211(5048):481–483.
- Beaupré SR, Druffel ERM. 2009. Constraining the propagation of bomb-radiocarbon through the dissolved organic carbon (DOC) pool in the northeast Pacific Ocean. *Deep Sea Research Part I: Oceanographic Research Papers* 56:1717–1726.
- Beaupré SR, Druffel ERM. 2012. Photochemical reactivity of ancient marine dissolved organic carbon. *Geophysical Research Letters* 39(18): L18602.
- Beaupré SR, Druffel ERM, Griffin S. 2007. A low-blank photochemical extraction system for concentration and isotopic analyses of marine dissolved organic carbon. *Limnology and Oceanography-Methods* 5(6):174–184.
- Broek TAB, Walker BD, Guilderson TP, McCarthy MD. 2017. Coupled ultrafiltration and solid phase extraction approach for the targeted study of semi-labile high molecular weight and refractory low molecular weight dissolved organic matter. *Marine Chemistry* 194:146–157.
- Druffel E, Griffin S, Walker BD. 2013. Total uncertainty of radiocarbon measurements of marine dissolved organic carbon and methodological recommendations. *Radiocarbon* 55(3):1135–1141.
- Griffin S, Beaupré SR, Druffel ERM. 2010. An alternate method of diluting dissolved organic carbon seawater samples for ^{14}C analysis. *Radiocarbon* 52(3):1224–1229.
- Lide DR. 2005. *CRC handbook of chemistry and physics*, internet version 2005. Boca Raton: CRC Press.
- Oppenländer T. 2007. Photochemical purification of water and air: advanced oxidation processes (AOPs)-principles, reaction mechanisms, reactor concepts. Weinheim: Wiley-VCH.
- Santos GM, Southon JR, Griffin S, Beaupré SR, Druffel ERM. 2007. Ultra small-mass AMS C-14 sample preparation and analyses at KCCAMS/UCI facility. *Nuclear Instruments & Methods in Physics Research B* 259(1):293–302.
- Walker BD, Druffel ERM, Kolasinski J, Roberts BJ, Xu X, Rosenheim BE. 2017. Stable and radiocarbon isotopic composition of dissolved organic matter in the Gulf of Mexico. *Geophysical Research Letters* 44(16):8424–8434.
- Walker BD, Griffin S, Druffel ERM. 2016. Effect of acidified versus frozen storage on marine dissolved organic carbon concentration and isotopic composition. *Radiocarbon* 59(3): 843–857.

- Walker BD, Xu X. 2019. An improved method for the sealed-tube zinc graphitization of microgram carbon samples and ^{14}C AMS measurement. *Nuclear Instruments & Methods in Physics Research B* 438:58–65.
- Williams PM. 1968. Stable carbon isotopes in the dissolved organic matter of the sea. *Nature* 219(5150):152.
- Williams PM, Oeschger H, Kinney P. 1969. Natural radiocarbon activity of dissolved organic carbon in north-east pacific ocean. *Nature* 224(5216):256.
- Xu X, Trumbore SE, Zheng S, Southon JR, McDuffee KE, Luttgen M, Liu JC. 2007. Modifying a sealed tube zinc reduction method for preparation of AMS graphite targets: reducing background and attaining high precision. *Nuclear Instruments & Methods in Physics Research B* 259(1):320–329.
- Xue Y, Ge T, Wang X. 2015. An effective method of UV-oxidation of dissolved organic carbon in natural waters for radiocarbon analysis by accelerator mass spectrometry. *Journal of Ocean University of China* 14(6):989–993.
- Zigah PK, McNichol AP, Xu L, Johnson C, Santinelli C, Karl DM, Repeta DJ. 2017. Allochthonous sources and dynamic cycling of ocean dissolved organic carbon revealed by carbon isotopes. *Geophysical Research Letters* 44(5):2407–2415.