

Benefits of Manometer in Non-Invasive Ventilatory Support

Rodrigo Silva Lacerda;¹ Fernando Cesar Anastácio de Lima;¹ Leonardo Pereira Bastos;² Anderson Fardin Vinco, MSc;² Felipe Britto Azevedo Schneider, MSc;² Yves Luduvico Coelho, MSc;² Heitor Gomes Costa Fernandes;² João Marcus Ramos Bacalhau;² Igor Matheus Simonelli Bermudes;² Claudinei Ferreira da Silva;¹ Luiza Paterlini da Silva;² Rogério Pezato, MD PhD^{1,3}

1. Department of Emergency Medical Services-Fire Academy, Escola Superior de Bombeiros-Polícia Militar do Estado de São Paulo, Franco da Rocha, Brazil
2. Inovelab, Laboratory of Engineering, Vitória, ES, Brazil
3. ENT Research Laboratory, Department of Otolaryngology-Head and Neck Surgery, Federal University of São Paulo, São Paulo, SP, Brazil

Correspondence:

Rogério Pezato, MD, PhD
Rua dos Otonis, 700 Piso Superior
Vila Clementino, 04025-002 São Paulo, SP,
Brazil
E-mail: pezatobau@ig.com.br

Conflicts of interest/funding: This manuscript had technical support from Vent-Logos (Vitória, ES, Brazil), the company responsible for the manufacture of automatic ventilators. The authors declare no conflicts of interest.

Keywords: airway; cardiac arrest; ventilatory support

Abbreviations:

AV: automatic ventilator
BVM: bag-valve-mask
CG: control group
CPR: cardiopulmonary resuscitation
ESB: Escola Superior de Bombeiros
MS: manometer support

Received: October 8, 2016

Revised: December 27, 2016

Accepted: January 28, 2017

Online publication: July 26, 2017

doi:10.1017/S1049023X17006719

Abstract

Introduction: Effective ventilation during cardiopulmonary resuscitation (CPR) is essential to reduce morbidity and mortality rates in cardiac arrest. Hyperventilation during CPR reduces the efficiency of compressions and coronary perfusion.

Problem: How could ventilation in CPR be optimized? The objective of this study was to evaluate non-invasive ventilator support using different devices.

Methods: The study compares the regularity and intensity of non-invasive ventilation during simulated, conventional CPR and ventilatory support using three distinct ventilation devices: a standard manual resuscitator, with and without airway pressure manometer, and an automatic transport ventilator. Student's t-test was used to evaluate statistical differences between groups. P values < .05 were regarded as significant.

Results: Peak inspiratory pressure during ventilatory support and CPR was significantly increased in the group with manual resuscitator without manometer when compared with the manual resuscitator with manometer support (MS) group or automatic ventilator (AV) group.

Conclusion: The study recommends for ventilatory support the use of a manual resuscitator equipped with MS or AVs, due to the risk of reduction in coronary perfusion pressure and iatrogenic thoracic injury during hyperventilation found using manual resuscitator without manometer.

Lacerda RS, de Lima FCA, Bastos LP, Vinco AF, Schneider FBA, Coelho YL, Fernandes HGC, Bacalhau JMR, Bermudes IMS, da Silva CF, da Silva LP, Pezato R. **Benefits of manometer in non-invasive ventilatory support.** *Prehosp Disaster Med.* 2017;32(6):615-620.

Introduction

Cardiopulmonary resuscitation (CPR) maneuvers are essential in the management of cardiac arrest. To increase the odds of early return of spontaneous circulation and, consequently, reduce downtime and mortality rate, CPR technique should be as effective as possible.

According to US resuscitation guidelines, one important aspect of CPR is that excessive ventilation should be avoided,¹ as hyperventilation is detrimental to patient survival rates.²

Hyperventilation during CPR is believed to reduce the efficiency of compressions and hinder venous return, as well as prevent the generation of negative intrathoracic pressure.²

In preclinical studies with a porcine model, the ratio of surviving animals declined to one out of seven from six out of seven when respiratory rate was increased to 30 versus 12 breaths/minute.³ This increase in ventilatory rate leads to a reduction in coronary perfusion pressure (regarded as the best hemodynamic marker of CPR efficacy), consequently reducing the odds of survival.³

To ensure that effective ventilation is achieved during CPR requires understanding which factors limit the proper performance of CPR technique by responders.

The present study was developed to compare the regularity and intensity of non-invasive ventilation during simulated, conventional CPR and ventilatory support performed by trained Basic Life Support responders using three distinct ventilation devices: a standard manual resuscitator (bag-valve-mask, BVM); a BVM equipped with a digital airway pressure manometer; or an automatic transport ventilator.




Figure 1. Manikin Configuration.

Methods and Materials

This study was carried out at Escola Superior de Bombeiros (ESB, the fire academy for the state of São Paulo, Brazil), located in the municipality of Franco da Rocha. The sample consisted of firefighters trained in prehospital emergency care.

Due to the risk of cardiothoracic injury when testing resuscitation techniques in human subjects, this study was performed on manikins; no victims or patients were involved.

The study protocol was evaluated and approved by the Ethics Committee of the ESB Department of Emergency Medical Services.

Methods and Measurements

Three groups were established for comparison purposes: the control group (CG), in which responders used BVMs alone, both for ventilatory support and for CPR; the manometer support (MS) group, in which responders used BVMs equipped with digital manometers, both for ventilatory support and for CPR; and the automatic ventilator (AV) group, in which responders used automatic transport ventilators, both for ventilatory support and for CPR. After sample size calculation, 11 responders (nine men, two women) were randomly selected from a cohort of firefighters squad (trained Basic Life Support responders) to perform ventilations in the three different study groups.

The CPR protocol used in the study consisted of 30 compressions followed by two breaths (compression-to-ventilation ratio 30:2), whereas the ventilatory assistance protocol consisted of 12 ventilations/minute.

Three CPR training manikins (QCPR, Laerdal; Wappinger, New York USA) equipped with data-logging airway monitors were used to measure the following variables: pressure, flow, volume, ventilation rate, inspiratory time, expiratory time, and inspiratory/expiratory (I:E) ratio (Figure 1). Pressure and flow measurements for each of the three monitor modules were calibrated and adjusted against an airway monitor previously calibrated at a Brazilian Calibration Network laboratory accredited by the National Institute of Metrology, Quality, and Technology (Inmetro; Rio de Janeiro, Brazil).

The devices used in the CG and MS groups were Oxigel (São Paulo, Brazil) brand, adult-size manual resuscitators, consisting of a reusable, autoclavable, silicone reservoir (capacity 1800 ml) and a one-way pop-off valve (pressure limit 60 cmH₂O). The devices used for the AV group were time-cycled,




Figure 2. Mean Peak Inspiratory Pressure with Ventilatory Support.

Abbreviations: AV, automatic ventilator; CG, control group; MS, manometer support.

pressure-limited, automatic transport ventilators (VentPlus; Vitória, Brazil) with a fixed I:E ratio.

Both the manual resuscitators and automatic ventilators were connected to Ventcare (Buri, Brazil) face masks (with silicone-cushioned rim, inflation valve, and hook ring) secured with Vent 11 model silicone straps. The head strap was used to improve mask seal.

The manometer developed for mechanical ventilation monitoring in this study is a low-pressure, low-power, digital gauge that measures peak inspiratory pressure and minimum airway pressure.

Outcomes

To evaluate the regularity and intensity of non-invasive ventilation during conventional CPR and ventilatory support, it was assessed the peak inspiratory pressure and ventilatory rate.

Statistics

After collection of quantitative data from all respiratory cycles, means and standard deviations were calculated for each group. The generated data were analyzed in PASW Statistics Version 18 (SPSS Inc.; Chicago, Illinois USA). Student's *t*-test was used to evaluate statistical differences between groups. *P* values < .05 were regarded as significant.

Results

The average age of responders was 24 years old (ranging from 19 to 27).

Peak Inspiratory Pressure during Ventilatory Support

Comparison between the CG and MS group revealed a statistically significant reduction ($P = .006$) in peak inspiratory pressure in the MS group during ventilatory support alone (Figure 2; Table 1).

Comparison between the MS and AV groups revealed a significantly lower mean inspiratory pressure ($P = .043$) in the latter. As expected, comparison between the CG and AV group revealed a significantly higher ($P = .0005$) peak inspiratory pressure in the CG, in which airway pressure during ventilation was neither monitored nor limited.




Figure 3. Mean Ventilation Rate with Ventilatory Support. Abbreviations: AV, automatic ventilator; CG, control group; MS, manometer support.

Ventilatory Rate during Ventilatory Support

Ventilatory rate was significantly higher in the CG than in the MS ($P=.002$) or AV ($P<.001$) groups (Figure 3; Table 2). Comparison between the MS and AV groups revealed a lower ventilatory rate in the latter ($P=.041$).

Peak Inspiratory Pressure during CPR

Peak inspiratory pressure during CPR was significantly higher in the CG than in the MS or AV groups ($P<.001$ for both comparisons; Figure 4; Table 3). Comparison of peak inspiratory pressure during CPR between the MS and AV groups revealed higher values in the former ($P=.041$).

Discussion

During ventilatory support, both hyperventilation and excessive inspiratory pressure can increase the risk of organ injury and reduce the patient's odds of recovery.^{2,3}

In mammals, a tidal volume of $6.3 \text{ mL}\cdot\text{kg}^{-1}$ is considered appropriate.⁴ The main causes of inadequate artificial ventilation are excessively high inspiratory pressure and excessively high respiratory rate.

High inspiratory pressures can lead to pulmonary hyperinflation, hemodynamic instability, gastric insufflation, and increased risk of barotrauma and, consequently, pneumothorax.^{2,5} Gastric insufflation is reported in up to 50% of patients receiving non-invasive pressure support ventilation. The lower esophageal sphincter pressure, or gastroesophageal junction pressure, which prevents reflux of gastric contents into the pharynx and




Figure 4. Mean Peak Inspiratory Pressure with CPR. Abbreviations: AV, automatic ventilator; CG, control group; CPR, cardiopulmonary resuscitation; MS, manometer support.

insufflation of air into the gastrointestinal tract during ventilation, is approximately $20\text{--}25 \text{ cmH}_2\text{O}$. The odds of gastric insufflation during BVM ventilation increase five-fold when peak airway pressure increases from 20 to $30 \text{ cmH}_2\text{O}$.⁵ High ventilatory rates, in turn, reduce coronary perfusion pressure.^{3,5}

In the present study, wide variation and higher ventilatory rates were observed when non-invasive ventilation was provided using a standard manual resuscitator alone, both during CPR and during ventilatory support, than when an AV was used or when the manual resuscitator was equipped with an airway pressure manometer. In this group, inspiratory pressure was also higher than in the MS and AV scenarios, which suggests that ventilation with a BVM alone is associated with increased risk of complications compared to ventilation with a manometer-equipped manual resuscitator or, especially, with an AV.

Limitations

The limitations of this study were the use of manikins instead of humans and a small number of responders assessed; more accurate studies will be necessary to confirm the manuscript's results.

Conclusion

The study suggests that, in prehospital emergency care settings, ventilatory support should preferably be provided using manual resuscitators equipped with airway pressure manometers or better, using portable AVs, due to the risk of reduction in coronary perfusion pressure and iatrogenic thoracic injury during hyperventilation found using manual resuscitators without manometers.

References

- Dabrowska A, Telec W. [New guidelines of basic and advanced cardiopulmonary resuscitation and emergency cardiovascular care (ECC) American Heart Association (AHA)]. *Wiad Lek.* 2011;64(2):127-131.
- O'Neill JF, Deakin CD. Do we hyperventilate cardiac arrest patients? *Resuscitation.* 2007;73:82-85.
- Aufderheide TP, Sigurdsson G, Pirralo RG, et al. Hyperventilation-induced hypotension during cardiopulmonary resuscitation. *Circulation.* 2004;109(16):1960-1965.
- Schultz MJ, Haitsma JJ, Slutsky AS, Gajic O. What tidal volumes should be used in patients without acute lung injury? *Anesthesiology.* 2007;106:1226-1231.
- Ho-Tai LM, Devitt JH, Noel AG, O'Donnell MP. Gas leak and gastric insufflation during controlled ventilation: face mask versus laryngeal mask airway. *Can J Anaesth.* 1998;45(3):206-211.

		Mean Peak Inspiratory Pressure (cmH ₂ O)											
		CG				MS				AV			
Round/ Simulator	Group	Mean	SD	Samples	CI	Mean	SD	Samples	CI	Mean	SD	Samples	CI
1/2	2	30.89	2.49	81	0.46	24.55	1.38	24	0.48	22.89	0.2	22	0.07
1/3	3	34.51	3.26	66	0.67	18.84	2.17	29	0.69	20.49	0.31	27	0.10
2/1	4	28.73	1.26	52	0.29	19.25	1.16	37	0.32	20.68	0.2	24	0.07
2/2	5	28.14	1.78	73	0.35	25.89	1.18	18	0.48	22.97	0.3	26	0.10
2/3	6	18.00	0.82	39	0.22	21.33	1.13	29	0.36	17.76	3.9	37	1.08
3/1	7	27.58	1.40	72	0.27	20.76	1.16	30	0.36	20.19	0.21	28	0.07
3/2	8	36.36	2.35	62	0.50	26.31	1.24	52	0.29	22.54	0.8	31	0.24
3/3	9	26.28	1.58	49	0.38	23.40	0.59	67	0.12	22.75	1.49	36	0.42
4/1	10	23.11	1.72	57	0.38	18.04	1.40	26	0.47	17.42	2.34	32	0.70
4/2	11	29.54	1.89	61	0.40	24.76	1.44	26	0.48	23.37	0.41	22	0.15
4/3	12	24.19	1.8	47	0.44	24.26	1.8	26	0.60	22.66	0.35	45	0.09
OVERALL		27.94	5.15	11	2.81	22.49	2.96	11	1.62	21.25	2.12	11	1.16

Lacerda © 2017 Prehospital and Disaster Medicine

Table 1. Analysis of Peak Inspiratory Pressure Data from Different Ventilatory Support Scenarios
Abbreviations: AV, automatic ventilator; CG, control group; MS, manometer support.

		Mean Ventilation Rate (breaths/min)											
		CG				MS				AV			
Round/ Simulator	Group	Mean	SD	Samples	CI	Mean	SD	Samples	CI	Mean	SD	Samples	CI
1/2	2	38.35	2.57	81	0.48	11.59	1.20	24	0.42	9.18	0.08	22	0.03
1/3	3	31.56	2.48	66	0.51	14.92	1.23	29	0.39	11.23	0.17	27	0.06
2/1	4	25.46	1.55	52	0.36	18.39	1.8	37	0.50	10.31	0.11	24	0.04
2/2	5	35.39	2.39	73	0.47	9.38	1.36	18	0.56	11.61	0.28	26	0.09
2/3	6	19.18	1.81	39	0.49	15.15	3.17	29	1.00	15.39	1.31	37	0.36
3/1	7	35.85	1.59	72	0.31	15.79	3.91	30	1.21	11.57	0.18	28	0.06
3/2	8	30.86	1.88	62	0.40	25.97	2.10	52	0.49	11.98	0.6	31	0.18
3/3	9	24.46	2.41	49	0.58	33.31	1.41	67	0.29	14.83	0.18	36	0.05
4/1	10	28.77	2.50	57	0.55	12.56	0.85	26	0.28	10.62	0.07	32	0.02
4/2	11	30.65	1.73	61	0.37	12.68	1.00	26	0.33	9.29	0.16	22	0.06
4/3	12	26.24	1.81	47	0.44	13.10	0.97	26	0.32	15.74	0.13	45	0.03
OVERALL		29.71	5.65	11	3.09	16.62	7.05	11	3.85	11.98	2.33	11	1.28

Lacerda © 2017 Prehospital and Disaster Medicine

Table 2. Analysis of Ventilation Rate Data from Different Ventilatory Support Scenarios
Abbreviations: AV, automatic ventilator; CG, control group; MS, manometer support.

		Mean Peak Inspiratory Pressure (cmH ₂ O)											
		CG				MS				AV			
Round/ Simulator	Group	Mean	SD	Samples	CI	Mean	SD	Samples	CI	Mean	SD	Samples	CI
1/2	2	32.13	3.77	12	1.95	18.18	1.98	10	1.15	19.05	0.79	12	0.41
1/3	3	32.45	2.83	12	1.47	24.7	0.91	12	0.47	20.77	1.07	12	0.55
2/1	4	24.25	6.18	10	3.58	16.54	3.49	8	2.34	16.63	3.7	12	1.92
2/2	5	29.18	1.23	12	0.64	26.06	1.09	10	0.63	21.39	0.68	10	0.39
2/3	6	23.77	0.76	10	0.44	22.18	1.45	12	0.75	16.59	1.86	10	1.08
3/1	7	27.03	1.7	10	0.99	19.86	1.03	10	0.60	20.21	0.65	12	0.34
3/2	8	38.38	1.85	10	1.07	26.82	2.46	10	1.43	21.42	0.54	12	0.28
3/3	9	23.07	2.81	10	1.63	21.18	1.09	10	0.63	20.38	2.03	12	1.05
4/1	10	30.48	4.84	10	2.81	17.99	1.96	10	1.14	15.41	3.42	10	1.98
4/2	11	35.76	5.54	14	2.62	25.59	1.35	12	0.70	21.66	0.9	12	0.47
4/3	12	27.25	3.42	10	1.98	19.22	3.76	10	2.18	22.27	1.43	12	0.74
OVERALL		29.43	4.98	11	2.72	21.67	3.64	11	1.99	19.62	2.37	11	1.29

Lacerda © 2017 Prehospital and Disaster Medicine

Table 3. Analysis of Peak Inspiratory Pressure Data from Different CPR Scenarios
Abbreviations: AV, automatic ventilator; CG, control group; CPR, cardiopulmonary resuscitation; MS, manometer support.